EN

OPERATING INSTRUCTIONS

- Twenty to thirty minutes before use, remove the product box from the refrigerator, remove the blister pack from the box and allow the syringe to come to room temperature. Do not use if the sterile package appears to be damaged.
- Peel off the blister backing.** (Fig. 1)
- While holding the blister open side down, bend the blister and allow the syringe to fall gently onto the sterile surface. (Fig. 2)
Alternatively, hold the blister open side up and bend back the blister until the barrel's luer end is exposed. Gripping the luer end of the barrel, remove the syringe from the blister. **Do not remove the syringe from the plunger end.**
- Remove the tip cap** (Fig. 3.1, Fig. 3.2) **from the syringe and attach the cannula** (Fig. 4.1, Fig. 4.2). **Attention:** Do not apply pressure to the plunger rod while the cannula is being affixed. Verify that the cannula is properly locked to the Luer Lock Adaptor (LLA). Do not overtighten the LLA; this can lead to loosening of the LLA from the barrel.
- Apply gentle pressure to the plunger** in order to expel air from the cannula and verify that the syringe is operating properly. (Fig. 5)
- The syringe is ready for use.
- Do not reuse or resterilize.

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REP

Cannula 27G

Sterimedix Ltd.
Thornhill Road
North Moons Moat
Redditch
B98 9ND, UK

0123

STERILE

EO

EC

REP

Bausch & Lomb GmbH
Brunsbütteler Damm 166-173,
Berlin, Germany 13581

STERILE

EO

Fig. 2

Fig. 3.1

Fig. 3.2

Fig. 4.1

Fig. 4.2

Fig. 5

6309501203

BiOLon®

BiOLon®

	Caution, see instructions for use		Batch code
	See instructions for use		Use by (YYYY-MM-DD : year-month-day)
	Do not reuse		Catalog number
	Protect from light		Medical device
	Do not use if the packaging has been opened or damaged		Keep dry
	Sterilized by ethylene oxide		Temperature limit
	Sterilized by aseptic techniques		Do not resterilize
	Manufacturer		Double sterile barrier system

BiOLon®

MANUFACTURED BY

Nordic Group B.V.
Siriusdreef 41,
2132 WT Hofdorp,
The Netherlands

DISTRIBUTED BY

NORDIC
PHARMA

6309501203 Issue date: xx/20xx

CE 2797

BIOLON®

Sodium hyaluronate, 10mg/ml

FORMULA	
Each ml of BIOLON® contains:	
Sodium hyaluronate	10 mg
Sodium chloride	8.5 mg
Disodium hydrogen phosphate dodecahydrate	0.56 mg
Sodium dihydrogen phosphate dihydrate	0.05 mg
Water for injections	q.s.

DESCRIPTION

Sodium hyaluronate is a physiological material found in both animal and human tissues. It is a polysaccharide consisting of a repeating disaccharide of N-acetylglucosamine and sodium glucuronate, linked by alternating β-1,3 and β-1,4 glycosidic bonds. Chemically identical in all species, hyaluronate can be found in the vitreous and aqueous humor of the eye, in synovial fluid, in skin and in the umbilical cord. BIOLON® is a sterile 1% solution of highly purified, high molecular weight viscoelastic sodium hyaluronate in phosphate-buffered saline. When introduced into the anterior segment of the eye during surgical procedures, BIOLON® maintains a deep anterior chamber throughout surgery, facilitates adhesions and helps to reneel the vitreous face and corneal endothelium and surrounding tissues. Its viscoelasticity also helps to reneel the vitreous face and discourage formation of a post-operative flat chamber in posterior segment surgery. BIOLON® serves as a surgical aid in separating, moving, and holding tissues. BIOLON® creates a clear field of vision for inspection of the retina and for photocoagulation during and after surgery. It does not interfere with epithelialization and normal wound healing. Any traces of BIOLON® left in the anterior segment of the eye after surgery dissipate through Schlemm's canal within a week.

INDICATION
Eye surgery, including intraocular lens insertion, intracapsular and extracapsular lens extraction, glaucoma surgery, corneal graft surgery for accidental trauma, retinal detachment and vitreal replacement procedures.

CONTRAINDICATIONS

When used as recommended, there are no known contraindications.

PRECAUTIONS

Precautions are limited to those normally associated with the surgical procedures being performed, including excessive amounts of BIOLON® into the anterior or posterior segments of the eye may increase intraocular pressure. Increased postoperative intraocular pressure may also be caused by a pre-existing glaucoma condition or by compromised outflow and by operative procedures and sequelae thereto, including enzymatic zonolysis, absence of an iridectomy, trauma to filtration structures, and by blood and lenticular remnants in the anterior chamber. Because these factors vary from case to case and are difficult to predict, the following precautions are recommended: Do not overfill the eye with BIOLON®. Avoid using large amounts of BIOLON® for posterior segment procedures in aphakic diabetic patients. At the end of the surgery, all remaining BIOLON® should be removed by irrigation or aspiration. Carefully monitor the intraocular pressure, especially during the immediate postoperative period. If a significant increase in pressure is observed, treat appropriately. Avoid the trapping of air bubbles behind BIOLON®. BIOLON® is a very highly purified substance extracted from bacterial cells. Since the presence of minute quantities of impurities (proteins, etc.) cannot be totally excluded, the physician should be aware of immunological, allergic and other potential risks associated with the injection of biological substances. On rare occasions, viscoelastic products containing sodium hyaluronate have been observed to become slightly opaque or to form a slight precipitate upon instillation into the eye. The clinical significance, if any, of this phenomenon is not known. The physician should, however, be aware of this possibility, and should it be observed, the cloudy or precipitated material should be removed by irrigation and/or aspiration. Avoid re-use of the cannula. If re-use is necessary, rinse the cannula thoroughly to remove all traces of residual materials.

DOSAGE

Dosage varies with type of surgery. Usually a dose of 0.2 to 0.6 ml is injected into the anterior segment of the eye. Greater amounts are used in the posterior segment.

EXAMPLES OF SURGICAL APPLICATIONS

Cataract Surgery: Before lens extraction, BIOLON® is introduced into the anterior chamber to protect the corneal endothelium and to maintain a deep anterior chamber. Additional amounts may also be introduced before implantation of an intraocular lens, and to coat the artificial lens and the surgical instruments.

Glaucoma Surgery: BIOLON® is instilled before the trabeculectomy in order to reconstitute the anterior chamber. At the close of surgery, further instillation of BIOLON® may be required to achieve good subconjunctival filtration and to prevent tissue adherences.

Corneal Transplant Surgery: In preparation of the cornea, the anterior chamber is filled with BIOLON®. The donor graft should be placed on the surface of the solution and sutured into position. Additional amounts BIOLON® may be injected to maintain a deep anterior chamber. BIOLON® may also be used in the anterior chamber of the donor eye to protect the corneal endothelial cells of the graft prior to preparation, and to protect the exposed endothelial layer of the donor button during preparation of the recipient eye.

Retinal and Vitreous Surgery: BIOLON® should be introduced slowly into the vitreous cavity. By carefully directing the injection, BIOLON® can be used to separate membranes from the retina for safe excision and release of traction. BIOLON® also helps in moving tissues into desired positions. For example, it is helpful gently pushing back a detached retina, in unrolling a retinal flap, or in holding the retina against the sclera during reattachment.

INCOMPATIBILITIES

Mixing of quaternary ammonium salts such as benzalkonium chloride with sodium hyaluronate solutions results in formation of a precipitate. BIOLON® should not be administered through a cannula previously used with medical solutions containing benzalkonium chloride.

ADVERSE REACTIONS

BIOLON® is well-tolerated in the human eye. Transient rises of post-operative intraocular pressure have been reported in some cases. A causal relationship has not been established between BIOLON® use and postoperative inflammatory reactions (iritis), corneal edema and corneal decompensation.

INTERACTIONS

None currently known.

PRESENTATIONS

BIOLON® is supplied in sterile disposable syringes containing 0.5 ml or 1.0 ml solution.

SHELF LIFE

3 years.

STORAGE INSTRUCTIONS

Store at 2°-8°C. May be kept at 25°C for up to six months. Do not freeze. Protect from light.

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