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




6309501305

⚠	Caution, see instructions for use	LOT	Batch code
📖	See instructions for use	📅	Use by (YYYY-MM-DD ; year-month-day)
🚫	Do not reuse	REF	Catalog number
💡	Protect from light	MD	Medical device
🚫	Do not use if the packaging has been opened or damaged	🧴	Keep dry
STERILE EO	Sterilized by ethylene oxide	🌡	Temperature limit
STERILE A	Sterilized by aseptic techniques	🚫	Do not resterilize
🏭	Manufacturer	🔒	Double sterile barrier system

MANUFACTURED BY
Nordic Group B.V.
Siriusdreef 41,
2132 WT Hofdorp,
The Netherlands

DISTRIBUTED BY
NORDIC
PHARMA

6309501205 Issue date: xx/20xx

Cannula 27G
Bausch & Lomb GmbH
Brunsbuteiler Damm 166-173,
Berlin, Germany 13581

EC REP
Stermedix Ltd.
Thornhill Road
North Moons Moat
Redditch
B98 9ND, UK

STERILE EO
0123

OPERATING INSTRUCTIONS

1. Twenty to thirty minutes before use, remove the product box from the refrigerator, remove the blister pack from the box and allow the syringe to come to room temperature.

Do not use if the sterile package appears to be damaged.

2. Peel off the blister backing. (Fig. 1)

3. While holding the blister open side down, bend the blister and allow the syringe to fall gently onto the sterile surface.

Alternatively, hold the blister open side up and bend back the blister until the barrel's luer end is exposed. Gripping the luer end of the barrel, remove the syringe from the blister. **Do not remove the syringe from the plunger end.** (Fig. 2)

4. **Remove the tip cap** (Fig. 3.1, Fig. 3.2) **from the syringe and attach the cannula** (Fig 4.1, Fig. 4.2). **Attention:** Do not apply pressure to the plunger rod while the cannula is being affixed. Verify that the cannula is properly locked to the Luer Lock Adaptor (LLA). Do not overtighten the LLA; this can lead to loosening of the LLA from the barrel.

5. **Apply gentle pressure to the plunger** in order to expel air from the cannula and verify that the syringe is operating properly. (Fig. 5)

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FORMULA

Each ml of BIOLON® Prime contains:

Sodium hyaluronate	12 mg
Sodium chloride	8.5 mg
Disodium hydrogen phosphate dodecahydrate	0.56 mg
Sodium dihydrogen phosphate dihydrate	0.05 mg
Water for injections	q.s

DESCRIPTION

BIOLON® Prime is a sterile 1.2% solution of highly purified, high molecular weight viscoelastic sodium hyaluronate in phosphate-buffered saline. Sodium hyaluronate is a physiological material found in both animal and human tissues. BIOLON® Prime is a very highly purified product extracted from bacterial cells. It is a polysaccharide consisting of a repeating disaccharide of N-acetylglucosamine and sodium glucuronate, linked by alternating β-1,3 and β-1,4 glycosidic bonds. Chemically identical in all species, hyaluronate can be found in the vitreous and aqueous humor of the eye, in synovial fluid, in skin and in the umbilical cord.

When introduced into the anterior segment of the eye during surgical procedures, BIOLON® Prime maintains a deep anterior chamber throughout surgery, facilitates surgery and reduces trauma to the corneal endothelium and surrounding tissues. Its viscoelasticity also helps to repel the vitreous face and discourage formation of a postoperative flat chamber.

Any traces of BIOLON® Prime left in the anterior segment of the eye after surgery dissipate through Schlemm's canal within a week.

INDICATION

Surgery in the anterior segment of the eye.

CONTRAINDICATIONS

When used as recommended, there are no known contraindications.

PRECAUTIONS

Precautions are limited to those normally associated with the surgical procedures being performed.

Instilling excessive amounts of BIOLON® Prime into the anterior or posterior segments of the eye may increase intraocular pressure. Increased postoperative intraocular pressure may also be caused by a pre existing glaucoma condition or by compromised outflow and by operative procedures and sequelae thereto, including enzymatic zonulysis, absence of an iridectomy, trauma to filtration structures, and by blood and lenticular remnants in the anterior chamber. Because these factors vary from case to case and are difficult to predict, the following precautions are recommended: Do not overfill the eye with BIOLON® Prime.

At the end of surgery, any remaining BIOLON® Prime should be removed by irrigation or aspiration.

Carefully monitor the intraocular pressure, especially during the immediate postoperative period. If a significant increase in pressure is observed, treat appropriately.

Avoid the trapping of air bubbles behind BIOLON® Prime.

Since the presence of minute quantities of impurities (proteins, etc.) cannot be totally excluded, the physician should be aware of immunological, allergic and other potential risks associated with the injection of biological substances.

On rare occasions, viscoelastic products containing sodium hyaluronate have been observed to become slightly opaque or to form a slight precipitate upon instillation into the eye. The clinical significance, if any, of this phenomenon is not known. The physician should, however, be aware of this possibility, and should it be observed, the cloudy or precipitated material should be removed by irrigation and/or aspiration.

Avoid re-use of the cannula. If re-use is necessary, rinse the cannula thoroughly to remove all traces of residual materials.

DOSAGE

Dosage varies with type of surgery. Usually a dose of 0.2 to 0.7 ml is injected into the anterior segment of the eye. Greater amounts are used in the posterior segment.

INCOMPATIBILITIES

Mixing of quaternary ammonium salts such as benzalkonium chloride with sodium hyaluronate solutions results in formation of a precipitate. BIOLON® Prime should not be administered through a cannula previously used with medical solutions containing benzalkonium chloride.

ADVERSE REACTIONS

BIOLON® Prime is well-tolerated in the human eye. Transient rises of postoperative intraocular pressure have been reported in some cases. A causal relationship has not been established between BIOLON® Prime use and postoperative inflammatory reactions (iritis), corneal edema and corneal decompensation.

INTERACTIONS

None currently known.

PRESENTATIONS

BIOLON® Prime is supplied in sterile disposable syringes containing 0.8 ml solution.

SHELF LIFE

3 years.

STORAGE INSTRUCTIONS

Store in a cold dark place (2°C-8°C). May be kept at 25°C for up to one month. Protect from freezing.