		RNC 2621
PRODUCT	LACRIFILL	
TERRITORY	UK FI SE DK NO CZ	
LANGUAGES	6 Languages EN FI SV DA NO CZ	
ART. CODE	LACN16B	
ARTICLE	Leaflet	
FORMAT	900 x 400 mm	
TEMPLATE	Tonnelier	
PRINTING INFORMATIONS		
COLORS		
1 COLOR BLACK		
FONT	Arial narrow	
MIN. FONT SIZE	9 pt	

Text update in line with the CA IFU
 V1 - 12/03/2026 (text updates / Section V)
 V2 - 20/03/2026 (text amendments)
 V3 - 24/03/2026 (text amendments)
 V4 - 30/03/2026 (date format)
 V5 - 31/03/2026 (SE date format change)

INSTRUCTIONS FOR USE

LACRIFILL® Canalicular Gel
 Disposable syringe
 Sterile Canalicular gel 20mg/ml

I. DESCRIPTION & COMPOSITION

Lacrifill® is sterile, colorless gel composed of non-animal hyaluronic acid (HA) cross-linked with 1,4-butanediol diglycidylether (BDDE) and reconstituted in a physiological buffer. The product is supplied in 1ml syringe (filling volume: 0.6ml), sterilized using moist heat sterilization process.

The pre-filled syringe is assembled with a plunger rod and a finger rest to allow gel placement. It is then packed in a secondary packaging (blister) to protect it from damage during transportation and storage and corresponds to one unit.

Lacrifill® is available in boxes of 1 or 10 units. Lacrifill® does not include any accessory. Qualitative and quantitative components of the injected gel:

Ingredients	Content per g
Crosslinked sodium hyaluronate	20 mg
Sodium chloride	9 mg
Potassium dihydrogen phosphate, KH ₂ PO ₄	0.02 mg
Dihydrogen hydrogen phosphate dihydrate, Na ₂ HPO ₄ ·2H ₂ O	0.22 mg
Water for Injection (WFI)	ad 1 g

II. INTENDED USE

a/ Intended purpose and indication

Lacrifill® is designed to temporarily block tear drainage by the occlusion of the lacrimal canaliculus of one or both eyelids in a given patient, thus maintaining lubricating tears on the surface of the eye.

Lacrifill® is intended for use, for up to 6 months, in adult patients experiencing dry eye symptoms such as redness, burning, reflex tearing, itching or foreign body sensations, which can be relieved by blocking of the canalicular system.

b/ Intended users

Lacrifill® shall only be administered by appropriately trained healthcare professionals who are qualified or accredited in accordance with national law, typically ophthalmologists or optometrists.

c/ Target population

Lacrifill® is used in adults experiencing dry eye symptoms such as redness, burning, reflex tearing, itching or foreign body sensations.

d/ Use environment

Lacrifill® is intended to be used in clinics/hospitals or private practices.

III. CONTRAINDICATIONS

Lacrifill® is contraindicated for patients with lacrimal outflow obstruction, inflammation of the ocular surface not associated with eye, active ocular or periorbital infection (for example dacryocystitis, canaliculitis, conjunctivitis, keratitis) and those who are allergic to hyaluronic acid.

IV. MODE OF ACTION

Lacrifill® is administered into the lower canaliculus of the patient's eyelid using a commercially available lacrimal cannula attached to the 0.6ml pre-filled Lacrifill® syringe. Lacrifill® is designed to temporarily block (for up to 6 months) tear drainage by the occlusion of the lacrimal canaliculus in a given patient, thus maintaining lubricating tears on the surface of the eye.

V. PREREQUISITES BEFORE AND OPERATING INSTRUCTIONS

FOR SINGLE PATIENT USE ONLY

Use of sterile gloves is recommended.

1. A sterile, disposable or reusable 25–27-gauge lacrimal cannula with Luer-lock connection may be used with Lacrifill® injection device. If a reusable lacrimal cannula is to be used, the cannula should be cleaned and sterilized prior to use for insertion of the gel consistent with the cannula's instructions for use and institutional procedures. Depending on the use of 25 or 27-gauge lacrimal cannula, the extrusion force will be slightly different, requiring more force to provide the same amount of gel using a 27-gauge lacrimal cannula.

2. Prior to punctal irrigation and insertion of Lacrifill®, a topical anesthetic agent should be instilled in the conjunctival sac for patient comfort. Punctal irrigation must demonstrate good lacrimal drainage prior to insertion of Lacrifill®.

3. To prepare for insertion of the gel, open the box, inspect the packaging for any signs of damage, remove the tray containing the gel-filled syringe, and peel off tray cover. Do not use if packaging has been compromised or device has been damaged.

4. Remove the cap from the gel-filled syringe and attach a sterile lacrimal cannula to the syringe (Figure 1).

II. INTENDED USE

a/ Intended purpose and indication

Lacrifill® is designed to temporarily block tear drainage by the occlusion of the lacrimal canaliculus of one or both eyelids in a given patient, thus maintaining lubricating tears on the surface of the eye.

Lacrifill® is intended for use, for up to 6 months, in adult patients experiencing dry eye symptoms such as redness, burning, reflex tearing, itching or foreign body sensations, which can be relieved by blocking of the canalicular system.

b/ Intended users

Lacrifill® shall only be administered by appropriately trained healthcare professionals who are qualified or accredited in accordance with national law, typically ophthalmologists or optometrists.

c/ Target population

Lacrifill® is used in adults experiencing dry eye symptoms such as redness, burning, reflex tearing, itching or foreign body sensations.

d/ Use environment

Lacrifill® is intended to be used in clinics/hospitals or private practices.


Figure 1: Attach cannula


Figure 2: Twist to lock cannula to syringe


Figure 3: Prime cannula with 0.1 mL gel


Figure 4: Place cannula tip in lower punctum


Figure 5: Deliver 0.2 mL gel into lower canaliculus

VI. WARNINGS

As with any procedure, use of Lacrifill® carries a risk of infection. Standard precautions associated with injectable materials should be followed. Should an infection occur, flush the gel from the canaliculus.

If the patient experiences irritation, infection, or epiphora after insertion of Lacrifill®, the gel should be removed.

In the pivotal clinical study, the safety and effectiveness of Lacrifill® was characterized over a 6-month period. The safety and effectiveness of the device for longer periods

VII. PRECAUTIONS

Lacrifill® is a sterile, colorless gel without visible particulate. If the contents of a syringe show signs of separation and/or appear cloudy or colored, do not use the syringe; notify Nordic Pharma immediately at lacrifill-vigilance@nordicpharma.com

Lacrifill® is provided STERILE for single patient use only for insertion in one or two lacrimal canaliculus during same session to avoid risk of infection.

Do not re-sterilize the syringe containing the gel. Do not insert more or less than 0.2 mL of gel into each lower punctum unless the cannular anatomy of a given patient requires it.

Use the device prior to the expiration date printed on the label.

Prior to use, inspect the package and device to verify that there is no damage to the device or package. Do not use the device if the device or packaging appears open or damaged.

Failure to comply with the cannula attachment instructions could result in cannula engagement and/or product leakage at the Luer Lock and cannula hub connection.

Prior to punctal irrigation and insertion of Lacrifill®, a topical anesthetic agent should be instilled in the conjunctival sac for patient comfort.

Do not insert more or less than 0.1 mL of gel before insertion in each lower canaliculus.

Inadequate or excessive suction force, or improper technique used to insert Lacrifill® into the lower canaliculus may cause patient discomfort.

After use, dispose of the syringe in accordance with accepted medical practice and applicable local, and state requirements.

If a reusable lacrimal cannula has been used, it should be cleaned and re-sterilized using standard procedures consistent with the cannula manufacturer's instructions for use and institutional procedures.

As with any punctal or cannular plug, Lacrifill® may enhance the effects of ocular medications used on the eye. Depending on the type of medications being used, dosage may need to be adjusted accordingly.

The clinical performance of Lacrifill® has not been established in the following conditions/patient populations:

- Pregnancy and breastfeeding.
- Pediatric patients.
- Use of ophthalmic cyclosporine (Restasis) within 6 months or lifitegrast (Xiidra) within 3 months.
- Corneal transplant.
- Ocular surgery (such as cataract surgery or LASIK) within six months.
- Current systemic infection, uncontrolled autoimmune disease, uncontrolled immunodeficiency disease.
- Significant corneal or conjunctival scarring, pterygium or nodular pinguecula; current ocular infection (except mild blepharitis), conjunctivitis or inflammation not associated with dry eye; anterior (epithelial) basement membrane corneal dystrophy or other clinically significant corneal dystrophy or degeneration; history of ocular herpetic infection; evidence of keratoconus; lid or lacrimal cancer.
- Active severe systemic allergy, seasonal allergies, rhinitis or sinusitis requiring treatment (i.e., antihistamines, decongestants, oral or aerosol steroids).
- Use of steroids, including administration by systemic or topical ocular routes.
- Ectropion.

VIII. POTENTIAL UNDESIRABLE SIDE EFFECTS

- (Eyelid) pain
- Allergic reactions
- Blepharitis

Canaliculitis

Conjunctivitis

Corneal ulceration

Dacryocystitis

Distal migration

Epiphora

Foreign Body Sensation

Itching

Local infection

Local irritation

Loss / Partial extrusion

Punctal granuloma

Swelling

IX. HOW LACRIFILL IS SUPPLIED AND STERILE CONDITIONS

Lacrifill® is supplied in a syringe containing 0.6 mL of gel for post-operative use on one side of a single patient during the same session. The contents of the syringe are sterile by moist heat method. Do not re-sterilize. Do not use if packaging is open or damaged.

X. EXPIRATION DATE, AND STORAGE, HANDLING AND DISPOSAL

Lacrifill® has a shelf life of 36 months. Lacrifill® must be used prior to the expiration date on the package. Store at a temperature of +5°C up to +25°C. Protect from sunlight. Refrigeration is not required.

After opening, discard any unused material of Lacrifill® according to the institutional policies.

XI. VIGILANCE AND POST-MARKET SURVEILLANCE

Any incident related to the device must be reported to Nordic Pharma. In the event of a serious incident that has occurred in relation to the device, it must be reported to Nordic Pharma and to the competent authority of the Member State in which the user and/or patient is established.

Notify Nordic Pharma immediately at lacrifill-vigilance@nordicpharma.com

Manufactured by: Nordic Pharma Nordic Group B.V. Sinusdreef 41, 2132 WT Hoofddorp, Netherlands

This IFU is available in electronic form on the legal manufacturer website: www.lacrifill.info

Latest approval date: 19-03-2026

of use has not been established. Therefore Lacrifill® should be used for no more than 6 months. After 6 months, Lacrifill® should be removed by thorough lacrimal irrigation.

When removing Lacrifill®, ensure that punctal irrigation is sufficiently thorough. If irrigation does not remove the gel, lacrimal probing or surgical exploration may be needed.

To minimize the risks of potential complications, Lacrifill® should only be inserted by health care practitioners who have appropriate training and experience, and who are knowledgeable about the anatomy at and around the punctum and canaliculus.

Obstruction of the canaliculus should be determined prior to insertion of Lacrifill®. Do not insert Lacrifill® if the canaliculus is obstructed. Care should be used not to perforate the canaliculus during the insertion of Lacrifill®. Perforation may cause pain and increased risk of infection. If perforation occurs, delay the insertion of the gel until the wound heals.

VII. PRECAUTIONS

Lacrifill® is a sterile, colorless gel without visible particulate. If the contents of a syringe show signs of separation and/or appear cloudy or colored, do not use the syringe; notify Nordic Pharma immediately at lacrifill-vigilance@nordicpharma.com

Lacrifill® is provided STERILE for single patient use only for insertion in one or two lacrimal canaliculus during same session to avoid risk of infection.

Do not re-sterilize the syringe containing the gel. Do not insert more or less than 0.2 mL of gel into each lower punctum unless the cannular anatomy of a given patient requires it.

Use the device prior to the expiration date printed on the label.

Prior to use, inspect the package and device to verify that there is no damage to the device or package. Do not use the device if the device or packaging appears open or damaged.

Failure to comply with the cannula attachment instructions could result in cannula engagement and/or product leakage at the Luer Lock and cannula hub connection.

Prior to punctal irrigation and insertion of Lacrifill®, a topical anesthetic agent should be instilled in the conjunctival sac for patient comfort.

Do not insert more or less than 0.1 mL of gel before insertion in each lower canaliculus.

Inadequate or excessive suction force, or improper technique used to insert Lacrifill® into the lower canaliculus may cause patient discomfort.

After use, dispose of the syringe in accordance with accepted medical practice and applicable local, and state requirements.

If a reusable lacrimal cannula has been used, it should be cleaned and re-sterilized using standard procedures consistent with the cannula manufacturer's instructions for use and institutional procedures.

As with any punctal or cannular plug, Lacrifill® may enhance the effects of ocular medications used on the eye. Depending on the type of medications being used, dosage may need to be adjusted accordingly.

The clinical performance of Lacrifill® has not been established in the following conditions/patient populations:

- Pregnancy and breastfeeding.
- Pediatric patients.
- Use of ophthalmic cyclosporine (Restasis) within 6 months or lifitegrast (Xiidra) within 3 months.
- Corneal transplant.
- Ocular surgery (such as cataract surgery or LASIK) within six months.
- Current systemic infection, uncontrolled autoimmune disease, uncontrolled immunodeficiency disease.
- Significant corneal or conjunctival scarring, pterygium or nodular pinguecula; current ocular infection (except mild blepharitis), conjunctivitis or inflammation not associated with dry eye; anterior (epithelial) basement membrane corneal dystrophy or other clinically significant corneal dystrophy or degeneration; history of ocular herpetic infection; evidence of keratoconus; lid or lacrimal cancer.
- Active severe systemic allergy, seasonal allergies, rhinitis or sinusitis requiring treatment (i.e., antihistamines, decongestants, oral or aerosol steroids).
- Use of steroids, including administration by systemic or topical ocular routes.
- Ectropion.

VIII. POTENTIAL UNDESIRABLE SIDE EFFECTS

- (Eyelid) pain
- Allergic reactions
- Blepharitis

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